

# C-SMART PILOT Endo motor

Large screen, but also wireless  
Simple, but also rich functional  
Split Endo motor & Apex locator

Three working mode



## Model: C-SMART

## Endo motors

### Feature:

- 1-Concise technology host, large 5-inch touch screen, rich information display, simple and easy to use.
- 2-The handle is wirelessly connected, without tedious connection lines.
- 3-The handle realizes the function measurement while expanding.
- 4- The contra angle can rotate freely at 300°, which can meet the treatment needs of any position.
- 5- Detachable integral insulation contra angle, no need insulation sheath.
- 6- The contra angle has built-in LED lighting for clearer vision.
- 7-The contra angle has a built-in root canal measurement connection, no need external snap.
- 8- Automatic calibration mode provides more accurate speed and torque.
- 9-Efficient and accurate motor drive and multi-frequency root canal measurement technology, can achieve root canal measurement and preparation at the same time
- 10- Built-in mainstream root canal expansion system, no need to set, just select and use
- 11-Imported special motor, quiet and high torque.
- 12- Wireless charging of the handle, the host and the handle can be charged at the same time.

### Technical:

- 1-Rotation Speed range: 150~1000rpm
- 2-Torque: 0.6~3.9N.cm
- 1-Adaptor working voltage: Input: AC100~240V, 50/60Hz  
Output: DC10V, 1.5A

Order name: **c-smart pilot**

This product has been applied layout-design patent, Counterfeit must be investigated.



The contra angle can rotate freely at 300°



Large 5-inch touch screen